



Metal Primer

A Metal Primer for in situ Repairs of Damaged
or Worn Veneers on Cast Restorations
Catalog No. 70-10

CONTENTS: 15cc Primer

INSTRUCTIONS FOR USE: Using a fine brush, apply a thin layer of *CuRAY-PriMet*® on washed, roughened, and blow-dried exposed metal surface. Wait one minute. Cure for ten seconds with visible light curing instrument. Follow with the application of light-cured restorative.

HELPFUL HINTS: The use of paint-on (flowable) light-cured restoratives, such as Scientific Pharmaceuticals' *CuRAY-Match*®, is strongly preferred over conventional filling restoratives.

Extend the restoration over a relatively large area of porcelain ground down to feather edge around the area to be restored.

STORAGE AND SHELF-LIFE: Store at temperatures not exceeding 75°F (24°C). When stored under such conditions, the material has a shelf-life of two years.

Quality Management System Certified to
ISO 9001 and ISO 13485
CE Marked Products

SCIENTIFIC PHARMACEUTICALS, INC.

3221 PRODUCER WAY • POMONA, CALIFORNIA USA 91768

PHONE: (800) 634-3047 • (909) 595-9922

FAX: (909) 595-0331 • E-MAIL: scipharm@msn.com

WEBSITE: www.scipharm.com

Trust these trademarks...

 **CuRAY-FIL**[®]

GLASION[®]

LUSTER[®]

 **CuRAY-MATCH**[®]

TEMPO-REX[®] **MICROJOIN**[®]

CEMPER[®] **DIAMOND LUSTER**[®]

 **CuRAY-SUPPORT**[®]

FloroPhos[®]

P CURE-ON-TOUCH **L**[®] **UNTIL**[®]

NEWCORE[™]

FLUROSEAL[®]

 **CuRAY - PriMet**[®]

FASTFIL[®]

SCI-SPAN[®]

LC SUPPORT[®]

 **CuRAY - II**[®]

SUPPORT[®]