



Manufactured in U.S.A. by:

SCI-PHARM®

Sci-Pharm Catalog Nos. 54-01/56-01
Light Cure Fissure Sealer and Bonding Agent

TECHNICAL BULLETIN Instructions

Light Cure Fissure Sealer and Bonding Agent

CONTENTS: Cat. No. 54-01 (Clear Kit)

- Cat. No. 54-011: Sealer (Clear), 6cc x 2 bottles
- Cat. No. 50-065: Enamel Conditioner, 13cc
- Mixing Wells
- Disposable brushes & Handle

CONTENTS: Cat. No. 56-01 (Opaque Kit)

- Cat. No. 56-011: Sealer (Opaque), 6cc x 2 bottles
- Cat. No. 50-065: Enamel Conditioner, 13cc
- Mixing Wells
- Disposable brushes & Handle

© Copyright 1985 Scientific Pharmaceuticals, Inc.

Cat. No. 54-01




OUTSTANDING FEATURES OF THE MATERIAL

- Versatility - Two versions available: The Clear (fluorescent) version is preferred for use as a Bonding Agent to improve adaption of composite restoratives to cavity margins and to provide retention in "adhesive" Class IV restorations; the Opaque version is preferred for use as a Pit and Fissure Sealer in a caries prevention program.
- Hydrophobicity - Low water absorption assures resistance to staining and durability in the oral environment.
- Ease in application - Simple dispensing system makes application easy and economical.
- No Refrigeration - Does not require refrigeration; shelf life is two years.

CAUTION

GENERAL INFORMATION

The use of fissure sealers for reducing the incidence of occlusal caries is well described in dental literature. The light-curing technique has made the application of sealers easier and more time efficient. Our Light-Cure Fissure Sealer and Bonding Agent represents an advanced formula that features low shrinkage, high stain resistance, low water sorption, outstanding depth of cure and long shelf-life with no need for refrigeration.



INSTRUCTIONS FOR USE

1. After prophylaxis, wash and dry the tooth (teeth) and apply Enamel Conditioner with a mini sponge, brush, cotton pledget, or other kind of an applicator on enamel only. **Caution:** Avoid contact with soft tissue or exposed dentin. If accidental spill occurs, wash immediately.
2. Wait one minute. Wash the tooth and evacuate. Dry with oil-free air or any commercial dental drying agent. Properly conditioned (etched) teeth should have a dull, chalky appearance. Highly mineralized teeth may require an additional two minute etch to achieve this effect.

If using the Opaque Sealer: *Before first use and before application if the material has not been used for more than one month, shake bottle well for approximately 15 seconds in order to ensure uniform dispersion of the opaquing pigment*

3. Dispense sealer into a mixing well. With a brush, apply a thin layer over the dry, etched enamel. Apply enough material to fill the fissure without, however, interfering with the occlusion. Cure with any commercial light for 20 seconds when used as a fissure sealer or for 10 seconds when used as a bonding agent.

Important: Maintain a dry field at all times during application and curing.

4. After the sealer sets, wipe the surface with a cotton pledget in order to remove the very thin film of uncured material remaining on the surface (the result of oxygen inhibition). Rinse and check the occlusion.

Federal law restricts this device to sale by or on the order of a dentist.

HELPFUL HINTS

1. The sealer is easily visible under UV-light when the clear (fluorescent) version is used. The opaque version permits easy detection of the sealer under any examination light.
2. Re-examination every six months after application of Fissure Sealer is recommended. In order to reapply, polish the tooth with a sanding disc and apply sealer following the procedure described above.
3. When sealing fissures, preferably work with one quadrant at a time.

STORAGE AND SHELF LIFE

Store at temperatures not exceeding 75°F (24°C). When stored under such conditions, the material has a shelf-life of 2 years.

For technical information, call or write:



Quality Management System Certified to

**ISO 9001 and ISO 13485
CE Marked Products**

SCIENTIFIC PHARMACEUTICALS, INC.

3221 PRODUCER WAY • POMONA, CALIFORNIA USA 91768

PHONE: (800) 634-3047 • (909) 595-9922

FAX: (909) 595-0331 • E-MAIL: scipharm@msn.com

WEBSITE: <http://www.scipharm.com>