



Manufactured in U.S.A. by:

**SCI-PHARM** 

**BAND SECURE**™

Light-Cured Orthodontic Band Cement  
Sci-Pharm Catalog No. 62-011

**TECHNICAL BULLETIN**  
Instructions

**BAND SECURE**™

**Orthodontic Band Cement with Fluoride**

**KIT CONTAINS:**

- Cement in Syringe-Type Dispenser, 5g
- Instructions

Cat. No. 62-011

**INDICATIONS FOR USE:**

For use as a direct bonding orthodontic adhesive for metal, plastic and ceramic brackets.

**CONTRAINDICATIONS:**

This product should not be used where patients have known hypersensitivity to methacrylate monomers.

**CAUTION: Enamel Conditioner contains Phosphoric Acid.**



**OUTSTANDING FEATURES  
OF THE MATERIAL**

- Virtually unlimited working time
- Ease of handling
- Economy of use – minimum waste
- Highly reliable
- Time savings, fast application
- Contains leachable fluoride
- Long shelf-life with no need for refrigeration

© Copyright 2007 Scientific Pharmaceuticals

**CAUTION**

**A. TOOTH AND BANDS PREPARATION**

- After removal of separators, scale gingival margins.
- Prophy teeth to be banded. Use prophy pastes that do not contain oils or petrolatums. Rinse and dry the teeth thoroughly.
- Trial fit bands. Clean bands before cementing. Use of bands with mechanically roughened or etched inside surfaces is recommended.

**B. CEMENTING**

- Isolate teeth to be banded.
- Squeeze the cement slowly from the syringe dispenser. Using a small spatula, or similar tool, cover the entire inside surface of the band with a thin layer of the cement. **While performing this procedure, shield the cement from direct exposure to intense light to avoid premature cure of the cement.** After dispensing, promptly replace the cap on the syringe.
- Make sure the tooth to be banded is dry and firmly seat the band. Remove excess cement with a scaler or cotton roll.
- Cure the cement from the occlusal side for 40 – 50 seconds (depending on light intensity) with conventional dental curing light, or 5 – 10 seconds when using high speed curing instrument.

**STORAGE AND SHELF LIFE**

Store at temperatures not exceeding 75°F (24°C). When stored under such conditions, the material has a shelf-life of two years.



*For technical information, call or write:*

*Quality Management System Certified to*  
**ISO 13485:2016**



**SCIENTIFIC PHARMACEUTICALS, INC.**

3221 PRODUCER WAY • POMONA, CALIFORNIA USA 91768

PHONE: (800) 634-3047 • (909) 595-9922

FAX: (909) 595-0331 • E-MAIL: [scipharm@msn.com](mailto:scipharm@msn.com)

WEBSITE: <http://www.scipharm.com>