



LUSTER®

Dental Polishing Paste For Composite Restoratives, Amalgams And Castings Cat. No. 50-05

KIT CONTAINS: 20g Paste (5 syringes of 4 grams each)

CAUTION: Federal Law restricts this device to sale by or on the order of a dentist.

OUTSTANDING FEATURES:

- Absolutely no drying or changes in consistency.
- Totally water soluble. Easily washes off teeth; cleans up and rinses off completely.
- Cohesive and non-splattering; will stay where you want it.
- Non-staining.
- Convenient and cost effective method of dispensing.

APPLICATION:

Finish in conventional manner using your favorite stones, discs or very fine diamonds, and then polish with Luster® using prophyl cups.

SHELF-LIFE AND STORAGE:

The shelf-life of Luster® is three years when stored at temperatures between 15°F (-10°C) and 83°F (30°C).

Quality Management System Certified to
ISO 9001 and ISO 13485
CE Marked Products

SCIENTIFIC PHARMACEUTICALS, INC.

3221 PRODUCER WAY • POMONA, CALIFORNIA USA 91768

PHONE: (800) 634-3047 • (909) 595-9922

FAX: (909) 595-0331 • E-MAIL: scipharm@msn.com

WEBSITE: www.scipharm.com

Trust these trademarks...

 **CURAY-FIL®**

GLASION®

LUSTER®

 **CURAY-MATCH®**

TEMPO-REX® MICROJOIN®

CEMPER® DIAMOND LUSTER®

 **CURAY-SUPPORT®**

FloroPhos®

PCURE-ON-TOUCH L® UNTIL®

NEWCORE™

FLUOROSEAL®

 **CURAY - PriMet®**

FASTFIL®

SCI-SPAN®

LC SUPPORT®

 **CURAY - II®**

SUPPORT®