

# DIAMOND LUSTER®

**Contents: 3g Polishing Paste**

**Diamond-filled Polishing Paste For Dental Ceramics  
(For Intra-oral & Extra-oral Use)**

**Cat. No. 55-05**

Scientific Pharmaceuticals' *Diamond Luster*® is an excellent material for restoring original smoothness to dental ceramics after making repairs, adjustments, etc. Designed for both intra- and extra-oral use, *Diamond Luster*® is based on high purity, non-toxic ingredients. A polyglycol base of controlled molecular weight prevents splattering and facilitates wash-out without staining the tissue. The abrasive consists of pure diamond powder, whose particle size distribution has been selected for best efficiency and optimum gloss.

The material is packaged in an economy size 3g syringe.

**METHOD OF APPLICATION:** The best way to apply *Diamond Luster*® is with felt polishing accessories.

**STORAGE AND SHELF-LIFE:** The shelf-life of *Diamond Luster*® is three years when stored at temperatures between 83°F (30°C) and 15°F (-10°C).

*Quality Management System Certified to*  
**ISO 9001 and ISO 13485**  
**CE Marked Products**

**SCIENTIFIC PHARMACEUTICALS, INC.**

3221 PRODUCER WAY • POMONA, CALIFORNIA USA 91768

PHONE: (800) 634-3047 • (909) 595-9922

FAX: (909) 595-0331 • E-MAIL: [scipharm@msn.com](mailto:scipharm@msn.com)

WEBSITE: [www.scipharm.com](http://www.scipharm.com)

**Trust these trademarks...**

 **CURAY-FIL®**

**GLASION®**

**LUSTER®**

 **CURAY-MATCH®**

**TEMPO-REX® MICROJOIN®**

**CEMPER® DIAMOND LUSTER®**

 **CURAY-SUPPORT®**

**FloroPhos®**

**PCURE-ON-TOUCH L® UNTIL®**

**NEWCORE™**

**FLUOROSEAL®**

 **CURAY - PriMet®**

**FASTFIL®**

**SCI-SPAN®**

**LC SUPPORT®**

 **CURAY - II®**

**SUPPORT®**